

Dip Slide Paddles



Intended Use

Hylabs products are agar slides specially designed to enable reliable, economic and timesaving monitoring of microbiological hygiene. Use directly on the work surface, dip into fluids or transfer from a conventional swab.

Product Specifications

Sensitivity:	Contact - 1 CFU/cm ² Dipping - 10 ³ CFU/cm ²
Testing time:	Total Count - 48 hours Coliformes and Generic <i>E. coli</i> - 18 hours Yeast and Mold - 48 hours
Storage:	Room temperature
Tests per kit:	20

Results



Ordering Information

Prod.#	Product description
9111	Hylabs Total Count
9114	Hylabs Yeast & Mold
9113	Hylabs Coliforms and Generic <i>E. coli</i>

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com