

NeoColumn®

For Aflatoxin DR (Direct Read)




Intended Use

NeoColumn Aflatoxin DR (Direct Read) is a immunoaffinity column designed for clean-up and concentration of sample prior to running on an HPLC system or direct read fluorescence for aflatoxin in commodities such as corn, corn by-products, peanuts and tree nuts.

Product Specifications

Lower limit of detection:	1 ppb
Testing time:	10 minutes
Recovery:	> 90% for B ₁ >80% for B ₂ , G ₁ , G ₂
Columns per kit:	50
Approvals:	USDA/GIPSA #2007-105 & #2007-105.1

Results



Ordering Information

Prod.#	Product description
8047	NeoColumn for Aflatoxin DR
8048	NeoColumn for Aflatoxin DR Developer - 50 mL
8049	NeoColumn for Aflatoxin DR Calibrator set

Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: foodsafety@neogen.com

Neogen Europe

+44 (0) 1292 525600 • E-mail: info@neogeneurope.com