

# Soleris®

## For Total Viable Count (TVC)



### Intended Use

Soleris Total Viable Count Vials are used as a screening method for detection of a variety of indicator and spoilage organisms.

### Product Specifications



Testing time: 4–24 hours  
Tests per kit: 100

### Materials Required But Not Provided

Available from Neogen

1. Soleris Automated System
2. Sterile PBS dilution bottles (BPB-99)
3. Adjustable pipette 1 mL

### Results



### Ordering Information

Prod.#	Product description
NB-100	Total Viable Count
NF-TVC	Non-Fermenting Total Viable Count

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)